

IHE Integration Statement			
Vendor	Product Name	Version	Date
Linkverse	Linkverse PACS: Linkverse Black-Storage Linkverse Black-Box Linkverse Black-Disc Linkverse Black-Vision	3	September 2009
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Implemented Integration Profiles	Actors Implemented	Options Implemented	
Scheduled Workflow (SWF)	Image Manager / Image Archive Image Display	None	
Patient Information Reconciliation (PIR)	Image Manager / Image Archive	None	
Consistent Presentation of Images (CPI)	Image Manager / Image Archive Image Display	None	
Portable Data for Imaging (PDI)	Portable Media Creator Portable Media Importer Image Display Display	Web content None None None	
Links to IHE			
In North America: www.rsna.org/IHE	In Europe: www.ihe-europe.net	In Japan: www.jira-net.or.jp/ihe-j	